**TONTARRA Medizintechnik GmbH**  
Surgical Instruments in Demand

### **HF-Chirurgie Zubehör**

HF-surgical accessories  
HF-chirurgical accessoires  
HF-quirúrgico accesorios  
HF-accessori per elettrochirurgia

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For more than 40 years TONTARRA has been recognized as a leader in quality and precision manufacturing of surgical instruments by continually implementing the latest advances in technology and production.

This success and the international recognition of the surgical skillfulness and accomplishment of TONTARRA Medizintechnik GmbH are based on our proven commitment to medical innovation, mastery in manufacturing and not least a strong focus on our continuing tradition of excellence.

The confidence in our medical programs and surgical lines as expressed for more than 40 years by our international customers and partners are a driving force to face and meet new challenges in surgical instrumentation.

Indispensable ‘tools’ for maturing new ideas into surgical instruments of highest precision include our state-of-the-art manufacturing facility and production processes and especially the ingenuity and strong commitment of our entire team.

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### Care instructions

Do not clean HF cables, adapters or electrode handles in an ultrasound bath. Do not sterilise HF cables, adapters or electrode handles and electrodes in hot air sterilisers. After use, wipe HF cable with an alcoholic surface disinfectant immediately. The HF cable or the handle may also be immersed in a combined cleaning and disinfecting solution; in this case, the service life may then be reduced due to oxidation of the contacts and crystallisation in the plugs. Observe the instructions of the chemical manufacturer and make sure that the agents used are compatible. Then rinse cable, connectors, pencil or electrode thoroughly with sterilised water. For concrete recommendations see separate instruction sheet which is included to the item or refer to catalogue page 86.

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







**advantage:**


- reusable
- current flow only partial
- only low HF-capacity: 30 till 60 Watt
- very low blistering
- to use in saline solution
- very good cut-/and vaporisation features
- suitable HF-cables

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3,0 m

5,0 m

286-215-30 international Models; 2-pins  
 286-215-50 international Models; 2-pins

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





Art.Nr.	Beschreibung   Description	Länge   Length	
	<b>bipolares Kabel für Aesculap</b> <b>bipolar cable for Aesculap</b> <b>câble bipolaire pour Aesculap</b> <b>cable bipolar para Aesculap</b> <b>cavo bipolare per Aesculap</b> <b>für Instrumentanschlüsse von:   for instrument plugs from:</b>		Umax: 0,8 KVp
287-500-50	Erbe ICC/ACC, Storz	5,0 m	
287-505-50	Berchtold; Martin	5,0 m	
287-510-50	US-Models; 2-pins	5,0 m	
287-515-50	international Models; 2-pins	5,0 m	
	<b>bipolares Kabel für Powerstar Ethicon</b> <b>bipolar cable for powerstar Ethicon</b> <b>câble bipolaire pour Powerstar Ethicon</b> <b>cable bipolar para Powerstar Ethicon</b> <b>cavo bipolare per Powerstar Ethicon</b> <b>für Instrumentanschlüsse von:   for instrument plugs from:</b>		Umax: 0,8 KVp
287-600-50	Erbe ICC/ACC, Storz	5,0 m	
287-605-50	Berchtold; Martin	5,0 m	
287-610-50	US-Models; 2-pins	5,0 m	
287-615-50	international Models; 2-pins	5,0 m	
	<b>bipolares Kabel für Storz Anschluss</b> <b>bipolar cable for Storz connection</b> <b>câble bipolaire pour fiche Storz</b> <b>cable bipolar para enchufe Storz</b> <b>cavo bipolare per attacco Storz</b> <b>für Instrumentanschlüsse von:   for instrument plugs from:</b>		Umax: 0,8 KVp
287-700-30	Erbe ICC/ACC, Storz	3,0 m	
287-700-50	Erbe ICC/ACC, Storz	5,0 m	
287-705-30	Berchtold; Martin	3,0 m	
287-705-50	Berchtold; Martin	5,0 m	
287-710-30	US-Models; 2-pins	3,0 m	
287-710-50	US-Models; 2-pins	5,0 m	
287-715-30	international Models; 2-pins	3,0 m	
287-715-50	international Models; 2-pins	5,0 m	

Art.Nr.	Beschreibung   Description	Länge   Length	
	<b>monopolares Kabel für Storz Resektoskop</b> <b>monopolar cable for Storz resectoscopes</b> <b>câble monopolaire pour Storz resectoscope</b> <b>cable monopolar para Storz resectoscope</b> <b>cavo monopolare per resettori Storz</b>		Umax: 3 KVp
287-800-30	monopolar cable 4 mm	3,0 m	
287-800-50	monopolar cable 4 mm	5,0 m	
	<b>Kabel für Neutralelektrode</b> <b>cable for return plates</b> <b>câble pour électrodes neutres</b> <b>cable para electrodos neutros</b> <b>cavo per piastre neutre monouso</b>		Umax: 3 KVp
287-805-30	Erbe; Martin	3,0 m	
287-805-50	Erbe; Martin	5,0 m	
	<b>Kabel für Einmal-Neutralelektrode</b> <b>cable for single use return plates</b> <b>câble pour électrodes neutres à usage unique</b> <b>cable para electrodos neutros desechable</b> <b>cavo per piastre neutre monouso</b>		Umax: 3 KVp
287-810-30	Erbe; Martin	3,0 m	
287-810-50	Erbe; Martin	5,0 m	
	<b>monopolares Kabel 4mm</b> <b>monopolar cable 4mm</b> <b>câble monopolaire 4mm</b> <b>cable monopolar 4mm</b> <b>cavo monopolare 4mm</b>		Umax: 3 KVp
287-815-30	Berchtold; Martin	3,0 m	
287-815-50	Berchtold; Martin	5,0 m	
	<b>monopolares Kabel 4mm</b> <b>monopolar cable 4mm</b> <b>câble monopolaire 4mm</b> <b>cable monopolar 4mm</b> <b>cavo monopolare 4mm</b>		Umax: 3 KVp
287-820-30	Erbe	3,0 m	
287-820-50	Erbe	5,0 m	
	<b>HF-Kabel nach Kundenwunsch, bzw. mit Sonderlängen können gerne angefertigt werden.</b> <b>Customers modifications are available upon request.</b>		

Art.Nr.	Beschreibung   Description	VPE   Unit
288-010-00	Bipolare Klemme; gebogen; gerieft bipolar clamp; curved; serrated pince bipolaire; courbe; striée pinza bipolar; curva; estriada pinza bipolare; curva; striata	1
	passende Kabel: 287-600-50 – 287-615-50 (Seite 78) suitable cables: 287-600-50 – 287-615-50 (page 78)	230 mm
	Bipolare Schere bipolar scissors ciseaux bipolaire tijera bipolar forbici bipolari	1
	gerade   straight droite   recta   retta	fein gebogen   curved courbe   curva   curva
288-020-18	288-021-18	180 mm
288-020-21	288-021-21	210 mm
288-020-23	288-021-23	230 mm
288-020-28	288-021-28	280 mm
	230 mm	
	Bipolares Kabel für bipolare Scheren bipolar cable for bipolar scissors câble bipolaire pour ciseaux bipolaire cable bipolar para tijera bipolar cavo bipolare per forbici bipolari	1
	für Instrumentanschlüsse von:   for instrument plugs from:	
288-030-30	US-Models; 2-pins	3,0 m
288-030-50	US-Models; 2-pins	5,0 m
288-035-30	Erbe ICC/ACC, Storz	3,0 m
288-035-50	Erbe ICC/ACC, Storz	5,0 m
288-040-30	Berchtold; Martin	3,0 m
288-040-50	Berchtold; Martin	5,0 m
288-045-30	international Models; 2-pins	3,0 m
288-045-50	international Models; 2-pins	5,0 m

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**Medical Device Directive (Europe)**  
**This device complies with the requirements of Directive 93/42 EEC concerning medical devices**

**1. Use**

HF cables connect electrode handles, electrodes and instruments for monopolar or bipolar coagulation to HF units. When used as intended, these products are more or less subject to wear, depending on the intensity of use. This kind of wear is caused by technical factors and cannot be avoided. Therefore carry out a thorough visual examination and a function test prior to each application of these products. In case a new product has externally visible defects or does not operate according to the instructions, contact the manufacturer or his responsible representative immediately.

**2. Operating**

Thoroughly examine the cable, plugs and electrode handle for irregularities on the surface. Do not use cables with a brittle or defective insulation. In case of doubt check functioning by means of a proper testing device.

**Connection to HF unit in monopolar mode**

First fix the neutral electrode ("return plate") to the patient, possibly on the upper arm or on the thigh. Ensure that the skin is free from hair and grease. Apply conductive gel evenly over the neutral electrode. The current paths inside the patient's body should be short and must proceed diagonally. Current paths must never run in transverse direction through the body or across the thorax. The patient must be insulated against all electrically conductive parts. Ground the operating table. Place the patient on a dry electrically insulating layer. Prevent skin areas touching each other. Place dry gauze in between. Then connect the neutral electrode, the foot switch and the "active" HF cable to the HF unit. In order to avoid interference with camera systems, do not guide HF cables in line with camera cables.

Never place HF cables directly on the

patient's skin. Never lay HF cables in loops. Then connect the HF cable to the instrument or to the electrode handle incl. the active electrode for coagulation and switch on the HF unit. The electrode with the proper shaft diameter (4mm resp. 2.4mm) has to be firmly attached to the handle tip.

Activation by yellow (finger) switch means: cutting current.

Activation by blue (finger) switch means: coagulating current.

Always refer to the instruction manual of the HF unit and to its general information about electrosurgical procedures. All monopolar products can be operated up to 8500 VP-P. (monopolar forceps up to 5500 VP-P).

**Connection to HF unit in bipolar mode**

In bipolar HF surgery the coagulation is restricted to the area between the two electrodes of the bipolar instrument. Consequently, a neutral electrode is not necessary. Besides, follow the recommendations regarding patient preparation in monopolar operating mode and the corresponding instructions of the HF unit. All bipolar products can be operated up to 2000 VP-P.

**3. Warnings**

**Power output of the HF unit:** Always select the lowest possible output of the HF unit for any procedure. If the coagulating capacity of the electrode is below normal, do not increase increase the output of the HF unit without a thorough preliminary check. Inspect for correct fit of the neutral electrode (in case of monopolar mode), correct attachment of all HF cables and plugs, correct activation of current by the finger- or foot-switches, undamaged insulation of the HF cables and instrument and check if the distal end of the electrode is clean and undamaged.

Before operating check if the signal transmission of the monitors works without interference (e.g. without noise); this check is done by pressing the yellow button (CUT) and the blue one (COAG).

Never place HF cables directly on the

**Endogenous burns**

Endogenous burns are those caused by high current density in the patient's tissue. Possible reasons are: The patient is inadvertently positioned in contact with electrically conductive parts. Direct contact between skin areas and HF cables may lead to electrical capacitance which in turn may cause burn.

**Exogenous burns**

Exogenous burns are those caused by the heat of ignited fluids or gases. They may also be caused by explosions. Possible reasons are among others: Ignition of skin cleaning agents and disinfectants, ignition of narcotic gases etc.

**Cardiac pacemakers**

Pacemakers may be damaged by electrosurgical current. Consult a cardiologist for the intervention. Never use electrosurgical current on outpatients fitted with pacemakers.

**4. Cleaning, disinfection and sterilisation**

**Caution:** Do not clean HF cables, adapters or electrode handles in an ultrasound bath. Do not sterilise HF cables, adapters or electrode handles and electrodes in hot air sterilisers.

After use, wipe HF cable with an alcoholic surface disinfectant immediately. The HF cable or the handle may also be immersed in a combined cleaning and disinfecting solution; in this case, the service life may then be reduced due to oxidation of the contacts and crystallisation in the plugs. Observe the instructions of the chemical manufacturer and make sure that the agents used are compatible. Then rinse cable, connectors, pencil or electrode thoroughly with sterilised water. For concrete recommendations see separate instruction sheet.

**5. Repairs**

**Caution:** Defective HF cables or electrode handles may not be repaired. Always replace a defective part with a new one.

• • •

**detachable bipolar coagulation forceps**

**Important Information - Please read before use !**

**Caution**

Please read all information contained in this insert. Incorrect handling and care as well as misuse can lead to premature wear of surgical instruments. This instrument is designed for use with the bipolar outputs of electrosurgical generators only. Do not exceed 60 Watts in the bipolar coagulation mode of generators.

**Description**

The Detachable Bipolar Coagulation Forceps are intended for use in general operative and laparoscopic surgical procedures. The device is intended to be passed through a 5mm laparoscopic cannula. Coagulation is achieved using electrosurgical energy under laparoscopic visualisation. This device is intended to be used with the outputs of compatible electrosurgical generators such as Erbe, Martin, Berchtold, Codman, Valleylab and similar generators. Do not exceed 60 Watts in the bipolar coagulation mode of generators.

**Indications for Use**

The Detachable Bipolar Coagulation Forceps are designed to manipulate, grasp and cut selected tissue. Bipolar Coagulation current may be selectively applied to the tissue as indicated. Suitable Cables:

**Note:** Bipolar coagulation instruments should be used only by individuals who are trained and licensed to use such devices.

Examples of such training and experience include:

Training through qualified residency program, surgical skills workshops, training programs offered by equipment manufacturers or preceptorship/surgical assistant training.

**Contraindications**

Not intended for contraceptive coagulation of the fallopian tube but may be used to achieve hemostasis following transection of the tube.

sparking in a flammable gas environment (i.e. explosive anesthetic gases). Organ perforation. Sudden massive hemorrhage.

**Contraindications to endoscopic procedures, not necessarily bipolar coagulation include:**

As identified in the Manual of Endoscopy available from the American Association of Gynecologic Laparoscopists. The presence of large pelvic or pelvic-abdominal masses, hypovolemic shock and severe cardiac decompensation. Also, intestinal obstruction and marked bowel distention, increase possibility of pelvic and abdominal adhesions.

A significantly elevated diaphragm contra-indicates the use of insufflation which may be necessary for proper surgical visualisation and may increase the chance of inadvertent bowel injury. Pelvic abscess, chronic pulmonary disease, diaphragmatic hernia, obesity, and septic peritonitis may exclude some patients from surgical consideration depending on severity of these conditions.

**Caution:** Please refer to the labeling and user manual for the electrosurgical generator for additional information on contraindications on electrosurgical or laparoscopic use.

Complications reported in the medical literature during laparoscopic surgical procedures include:

Hemorrhage, damage to surrounding soft tissue, leakage of bile or other secretions, infection (local and systemic), bowel perforation, damage to large blood vessels and/or neurological structures, inadvertently retained instruments and death.

**Adverse events reported while using bipolar electrosurgical devices include:**

Inadvertent activation with resultant tissue damage at the wrong site and/or equipment damage.

Fires involving surgical drapes and other combustible materials have been reported.

Alternate current pathways resulting in burns where the patient or physician or assistant is in contact with exposed metal.

Explosions caused by electrosurgical

Cautions & Warnings  
**US Federal Law restricts this device to sale by or on the order of a physician only.**

Special safety precautions should be observed when using electrosurgical instruments. Electrosurgical instruments can pose a significant shock, burn or explosion hazard if used improperly, incorrectly or carelessly.

Avoid touching or grounding electro-surgical instruments to non-insulated instruments, scopes, trocar sleeves, etc. All persons using such devices should be knowledgeable in the use and handling of laparoscopic instruments, coagulation equipment, their accessories and other related equipment.

Test all instruments, accessories and equipment prior to each use. Written Standard Operating Procedures for the cleaning, sterilisation, storage, inspection and maintenance of the instruments, accessories and equipment are recommended.

Do not use electrosurgical instruments on patients with pacemakers. Do not use in presence of flammable liquids or anaesthetics.

Electrosurgical generators used with these devices are designed to cause destruction of tissue and are inherently dangerous if operated improperly. Follow all safety precautions and instructions supplied by the manufacturer of the electrosurgical generator.

The electrode tip must always be in full view before activating power. Apply power only when electrode tip is in full contact with the tissue selected for coagulation.

Electrode tip must not come in contact with the laparoscope or other metal instruments during use.

Failure to observe these cautions and contraindications may result in injury, malfunction or other unanticipated

occurrences or events for the operator, staff and/or the patient.

#### Initial use of new instruments

Every instrument must be cleaned and sterilised before it is used for the first time.

The instrument was developed for sterilisation by autoclave and has shown good results using this method.

#### Inspection and functional check

It is very important to carefully examine each surgical instrument for breaks, cracks or malfunctions before use. It is especially essential to check areas such as blades, points, ends, stops and snaps as well as all movable parts.

Do not use damaged instruments. Never attempt to make repairs yourself. Service and repairs should be referred to trained qualified persons only. Refer questions about repair to the manufacturer or your biomedical engineering department.

#### Assembly and Operation

Assemble/disassemble forceps following the example assembly instructions and diagram.

Once correct assembled and connected to power, the device may be used in either the right or the left hand and any desired handle orientation to facilitate use and reduce hand discomfort and/or fatigue.

To close jaws: compress (grip) handle. To open jaws: gently release handle compression (grip).

Forceps jaws will not retract into the tube.

Never withdraw laparoscopic or surgical devices from the surgical site without direct observation or laparoscopic view to prevent inadvertent damage to adjacent tissue.

The Detachable Bipolar Coagulation Forceps are designed to provide concentrated grasping force or to cut at the forceps jaw and will cause some degree of tissue crushing where the jaws grasp or close. Tissue that is not to be clamped, crushed, cut or coagulated should not be grasped except with extreme care to avoid damage.

Coagulation forceps require the use of good surgical judgement and should not be used where crushing injury or coagulation is not acceptable. When retracting tissue with the forceps, the traction force applied to the tissue should be low enough to avoid tearing the selected tissue.

#### Cleaning and Maintenance

Every surgical instrument should be disinfected and thoroughly cleaned after each use. Proper cleaning, inspection and maintenance will help ensure correct function of the surgical instrument. Clean, inspect and test each instrument carefully. Sterilise all instruments before surgery. A good cleaning and maintenance procedure will extend the useful life of the instrument.

Special attention must be paid to slots, stops, ends, hollow tubes and other highly inaccessible areas. Check insulation, cables and connectors for cuts, voids, cracks, tears, abrasions, etc.

Do not use damaged instruments. This device can be disassembled. It must be cleaned in the disassembled state. (Note: Keep all parts together and protected from being misplaced, do not intermix with other or similar devices of other manufacture.)

Cleaning and rinsing must take place immediately after each use for best effect. Failure to clean promptly may result in adherent particles or dried secretions that may resist cleaning and complicate or resist future sterilisation. Instruments must be completely cleaned and rinsed of all foreign matter. Use warm water and a commercially available instrument pre-soak or cleaning agent. Enzymatic cleaners must be used to remove protein deposits. Follow the enzymatic cleaner's instructions; rinse thoroughly.

- Do not use corrosive cleaning agents (i.e. bleach). Cleaning solutions and rinses at or near a neutral pH (7.0) are best.
- Do not use abrasive cleaners.
- Only a soft bristle brush should be used.
- Immerse the entire device in deter-

gent and clean while soaking. Use a minimum of six strokes in and six strokes out with an instrument brush for all inside channels. Rinse with sterile deionized water. Repeat this procedure.

- Can be disinfected in the washing machine up to 203°F (95°C).
- Rinse thoroughly with distilled water.
- Prepare for storage and/or sterilisation.

After cleaning and rinsing, dry instruments completely and carefully with compressed air (including inside channels and highly inaccessible areas). Note: After cleaning and before sterilisation, treat all instruments with an oil which is considered as being physiologically safe, especially their blades, ends, stops, snaps and all movable parts.

#### Storage and Sterilisation

Instruments must be stored in a clean, dry, moisture free area. The instruments should be stored individually in their shipping carton or in a protective tray with partitions. Protects tips with cloth, gauze or tubing if stored in drawers.

Instruments and cables are reusable and meet AAMI standards for sterilisation.

We guarantee our products to withstand a minimum of twenty sterilisation cycles when sterilised according to the criteria listed below. Most devices will surpass this standard if a protocol is established for the proper care and inspection of medical devices at your facility.

Use steam autoclave sterilisation. Thoroughly clean instruments of all debris, tissue and foreign matter prior to sterilisation.

Follow the sterilizer manufacturer's instructions for operation and loading of steam autoclaves. There must be direct steam exposure to all surfaces of the instruments being sterilised including the internal surface and tubes channels.

Allow instrument to air cool to room temperature before use.

#### Standard sterilisation method

Use steam autoclave sterilisation only.

Standard autoclave cycle:

Steam sterilise at 270°F (132°C) for 15 minutes. Other time and steam temperature cycles may also be used. However, user must validate any deviation from the recommended time and temperature.

(Note: Contact the manufacturer of your steam autoclave to confirm appropriate temperatures and sterilisation times.)

Caution: Autoclave temperatures should not exceed 280°F (137°C); handles, insulation or other non-metallic parts may be damaged.

#### Handling

All surgical instruments must be handled with the greatest care when being transported, cleaned, treated, sterilised and stored.

This is especially true for blades, fine points and other sensitive areas. Surgical instruments corrode and their functions are impaired if they come into contact with aggressive materials. The instruments must not be exposed to acids or other aggressive cleaning agents.

#### Warranty

All products are guaranteed to be free from defects in material and workmanship at the time of shipping.

All of our products are designed and manufactured to meet the highest quality standards. We cannot accept any liability for failure of products which have been modified in any way from their originals, for any failure due to misuse or application which is not in accordance with the designer's intentions.

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Daimlerstraße 15

78573 Wurmlingen

Germany

Tel: +49(0) 74 61 / 9 65 76 -0

Fax: +49(0) 74 61 / 9 65 76 -26

[info@tontarra.de](mailto:info@tontarra.de)

[www.tontarra.de](http://www.tontarra.de)

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**TONTARRA MEDIZINTECHNIK GmbH**